

6. Table for Pt1000 probe

| Temperatura | R nominal (Ohm) |
|-------------|-----------------|
| -200        | 185,281         |
| -190        | 228,327         |
| -180        | 271,029         |
| -170        | 313,408         |
| -160        | 355,484         |
| -150        | 397,277         |
| -140        | 438,803         |
| -130        | 480,081         |
| -120        | 521,127         |
| -110        | 561,954         |
| -100        | 602,578         |
| -90         | 643,012         |
| -80         | 683,267         |
| -70         | 723,355         |
| -60         | 763,286         |
| -50         | 803,068         |
| -40         | 842,71          |
| -30         | 882,218         |
| -20         | 921,6           |
| -10         | 960,859         |
| 0           | 1000            |
| 10          | 1039,025        |
| 20          | 1077,936        |
| 30          | 1116,731        |
| 40          | 1155,411        |
| 50          | 1193,976        |
| 60          | 1232,426        |
| 70          | 1270,761        |
| 80          | 1308,981        |
| 90          | 1347,085        |
| 100         | 1385,075        |
| 110         | 1422,949        |
| 120         | 1460,709        |
| 130         | 1498,353        |
| 140         | 1535,882        |
| 150         | 1573,296        |

| Temperatura | R nominal (Ohm) |
|-------------|-----------------|
| 160         | 1610,595        |
| 170         | 1647,779        |
| 180         | 1684,848        |
| 190         | 1721,801        |
| 200         | 1758,64         |
| 210         | 1795,363        |
| 220         | 1831,972        |
| 230         | 1868,465        |
| 240         | 1904,843        |
| 250         | 1941,106        |
| 260         | 1977,254        |
| 270         | 2013,287        |
| 280         | 2049,205        |
| 290         | 2085,007        |
| 300         | 2120,695        |
| 310         | 2156,267        |
| 320         | 2191,725        |
| 330         | 2227,067        |
| 340         | 2262,294        |
| 350         | 2297,406        |
| 360         | 2332,403        |
| 370         | 2367,285        |
| 380         | 2402,052        |
| 390         | 2436,703        |
| 400         | 2471,24         |
| 410         | 2505,661        |
| 420         | 2539,968        |
| 430         | 2574,159        |
| 440         | 2608,235        |
| 450         | 2642,196        |
| 460         | 2676,042        |
| 470         | 2709,773        |
| 480         | 2743,389        |
| 490         | 2776,889        |
| 500         | 2810,275        |
| 510         | 2843,545        |

| Temperatura | R nominal (Ohm) |
|-------------|-----------------|
| 520         | 2876,701        |
| 530         | 2909,741        |
| 540         | 2942,666        |
| 550         | 2975,476        |
| 560         | 3008,171        |
| 570         | 3040,751        |
| 580         | 3073,216        |
| 590         | 3105,565        |
| 600         | 3137,8          |
| 610         | 3169,919        |
| 620         | 3201,924        |
| 630         | 3233,813        |
| 640         | 3265,587        |
| 650         | 3297,246        |
| 660         | 3328,79         |
| 670         | 3360,219        |
| 680         | 3391,533        |
| 690         | 3422,731        |
| 700         | 3453,815        |
| 710         | 3484,783        |
| 720         | 3515,637        |
| 730         | 3546,375        |
| 740         | 3576,998        |
| 750         | 3607,506        |
| 760         | 3637,899        |
| 770         | 3668,177        |
| 780         | 3698,34         |
| 790         | 3728,387        |
| 800         | 3758,32         |
| 810         | 3788,137        |
| 820         | 3817,84         |
| 830         | 3847,427        |
| 840         | 3876,899        |
| 850         | 3906,256        |